RADEX[™]



Direct-to-Digital General Radiographic System

Specially designed for ambulatory and outpatient general radiography exams.

The Hologic RADEX™ General Radiographic System is a direct-to-digital, fully integrated, imaging system for general purpose outpatient radiology. It utilizes a single digital detector to perform a wide variety of general radiography procedures, replacing both the wall and table Bucky in a conventional radiographic room. The system is designed for ambulatory patients providing high throughput, excellent image quality, and long-term reliability.

Built to optimize DirectRay direct-to-digital image technology

DirectRay technology directly captures and converts X-ray photons into a digital image, with no intermediate steps to compromise image quality.

The RADEX system consists of five major components:

- DirectRay® X-ray image capture system
- DirectRay Operator Console
- Wall-mounted U-arm
- X-ray generator
- Radiographic table





RADEX

Direct-to-Digital General Radiographic System

Product Specifications:

DirectRay® direct-to-digital X-ray image capture system

- · Digital detector array with full field, 14 x 17 inch/34 x 42 cm active image area
- 139 µm detector element pitch
- Amorphous selenium semiconductor X-ray absorber coating over thin-film transistor detector element matrix
- 3.6 cycles/mm detector Nyquist frequency
- Exposure factors available equivalent up to 400 speed screen-film systems
- Detail equivalent to highresolution screen-film systems
- Wide dynamic range–14-bit data captures; 12-bit output to PACS and film printers

DirectRay Operator Console

The Operator Console interfaces with the DirectRay Detector Array and the X-ray exposure equipment, while providing upstream and downstream connectivity to acquire and transmit patient, exam, and image data in digital format.

Wall-mounted U-arm

- · U-arm linked X-ray tube and
- Vertical travel: 43 inches/109
- Horizontal beam height, to floor: 20-63 inches, 51-160 cm
- Arm rotation: ±135°
- Motorized SID: 40-72 inches, 100-180 cm
- Digital SID display
- Electromechanical positioning brakes
- · Automatic collimation with manual override
- Bucky tilt: ±45°
- Bucky rotation for portrait and landscape array orientation
- Removable 10:1 grid
- 400,000 heat unit X-ray tube
- -0.6/1.2 mm focal spots
- 34/100 kW, 13.5° target
- Optional bracket provides 18" clearance from wall, allowing stretcher-based chest exams

Side View

Optional patient grip

X-ray generator system


- 80 kW high frequency output
 - 1000 mA @ 80 kVP
- 800 mA @ 100 kVP
- 500 mA @ 150 kVP
- -40 150 kVP
- 25 1000 mA
- Generator control integrated into DirectRay Operator Console
- Anatomically programmed radiography (APR) and automatic exposure control (AEC)
- Two-point and three-point

manual techniques and manual APR/AEC override capabilities

- Dual speed starter
- Line voltage range ± 10%
- 400 VAC, 480 VAC, 3 phase

Patient support device


- 300 lb/136 kg weight capacity
- Polycarbonate/plexiglass table
- Width: 24 inches/61 cm
- Length: 78 inches/199 cm
- Height: 30 inches/77 cm




Typical RADEX system room layout

Equipment legend					
Item	Description	Wt/lbs	Width	Depth	Height
001	RADEX - wall-mounted U-arm	900	Var.	40"	84"
002	RADEX - patient support device	110	78"	24"	Var.
003	X-ray generator cabinet	500	28"	43"	71"
004	DirectRay Operator Console	350	34"	20"	40"

Front View







Hologic Systems Division 300 Foster Street Littleton, MA 01460 USA Tel: 978.486.9681 Fax: 978.486.3426 www.hologic.com

Europe Horizon Park Leuvensesteenweg 510, BUS 31 1930 Zaventem, Belgium Tel: 32.2.711.4680 Fax: 32.2.725.2087

Asia, Pacific Rim and South America 35 Crosby Drive Bedford, MA 01730 USA Phone: 1.781.999.7300 Fax: 1.781.280.0668 www.hologic.com

The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications which must be met by Hologic, Inc. Characteristics of the products described in this publication may be changed at any time without notice.